

---

# **Blue Interface Documentation**

**Rachel Thomasson, Brent Yi**

**Oct 06, 2020**



---

## Contents

---

<b>Python Module Index</b>	<b>5</b>
<b>Index</b>	<b>7</b>



Blue Interface is a platform-agnostic Python API for controlling Blue robotic arms over a network connection.

It features:

- No dependency on ROS (or any particular version of Ubuntu)
- Easy connection to multiple robots
- Support for both Python 2 and 3
- Support for Mac, Windows, and Linux
- Support for Jupyter Notebooks

It's designed to be lightweight and easy-to-use! Sending a Blue “right” arm to its zero position, for example, is as simple as:

```
from blue_interface import BlueInterface

blue = BlueInterface(side="right", ip="127.0.0.1")
blue.set_joint_positions([0] * 7)
```

See [Github](#) for installation instructions and more usage examples.

---

**class** `blue_interface.BlueInterface` (*side, ip, port=9090*)

A Python interface for controlling the Blue robot through rosbridge.

#### Parameters

- **side** (*str*) – side of the arm, “left” or “right”
- **ip** (*str*) – The IP address of the robot, which by default should have a running rosbridge server.
- **port** (*int, optional*) – The websocket port number for rosbridge. Defaults to 9090.

**calibrate\_gripper** ()

Run the gripper position calibration process. This will automatically determine the gripper position by apply a closing torque and detecting when the gripper has fully closed.

**cancel\_gripper\_command** ()

Cancel current gripper command, halting gripper in current position.

**command\_gripper** (*position, effort, wait=False*)

Send a goal to gripper, and optionally wait for the goal to be reached.

#### Parameters

- **position** (*float64*) – gap size between gripper fingers in cm.
- **effort** (*float64*) – maximum effort the gripper with exert before stalling in N.

**disable\_control** ()

Set joint control mode to gravity compensation only.

**disable\_gripper** ()

Disables the gripper. The gripper will become compliant.

**enable\_gripper** ()

Enables the gripper. The gripper will begin to hold position.

**get\_cartesian\_pose** ()

Get the current cartesian pose of the end effector, with respect to the world frame.

**Returns** Pose in the form {"position": numpy.array([x,y,z]), "orientation": numpy.array([x,y,z,w]) defined with respect to the world frame.

**Return type** dict

**get\_gripper\_effort()**

Get the current effort exerted by the gripper.

**Returns** the gripper effort in N

**Return type** float64

**get\_gripper\_position()**

Get the current gap between gripper fingers.

**Returns** the gripper gap in cm.

**Return type** float64

**get\_joint\_positions()**

Get the current joint angles, in radians.

**Returns** An array of 7 angles, in radians, ordered from proximal to distal.

**Return type** numpy.ndarray

**get\_joint\_torques()**

Get the current joint torques.

**Returns** An array of 7 joint torques, in Nm, ordered from proximal to distal.

**Return type** numpy.ndarray

**get\_joint\_velocities()**

Get the current joint velocities.

**Returns** An array of 7 joint torques, in Nm, ordered from proximal to distal.

**Return type** numpy.ndarray

**gripper\_enabled()**

Check if gripper is enabled to take commands.

**Returns** True if enabled, False otherwise.

**Return type** bool

**inverse\_kinematics**(*position, orientation, seed\_joint\_positions=[]*)

Given a desired cartesian pose for the end effector, compute the necessary joint angles. Note that the system is underparameterized and there are an infinite number of possible solutions; this will only return a single possible one.

**Parameters**

- **position** (*iterable*) – A length-3 array containing a cartesian position (x,y,z), wrt the world frame.
- **orientation** (*iterable*) – A length-4 array containing a quaternion (x,y,z,w), wrt the world frame.
- **seed\_joint\_positions** (*iterable, optional*) – An array of 7 joint angles, to be used to initialize the IK solver.

**Returns** An array of 7 joint angles, or an empty array if no solution was found.

**Return type** numpy.ndarray

**set\_joint\_positions** (*joint\_positions*, *duration=0.0*, *soft\_position\_control=False*)

Move arm to specified position in joint space.

**Parameters**

- **joint\_positions** (*iterable*) – An array of 7 joint angles, in radians, ordered from proximal to distal.
- **duration** (*float*, *optional*) – Seconds to take to reach the target, interpolating in joint space. Defaults to 0.
- **soft\_position\_control** (*bool*, *optional*) – Use “software” position control, which runs position control loop at the ROS-level, rather than on the motor drivers. This should be rarely needed. Defaults to False.

**set\_joint\_torques** (*joint\_torques*)

Command joint torques to the arm.

**Parameters**

- **joint\_torques** (*iterable*) – An array of 7 joint torques, in Nm,
- **from proximal to distal.** (*ordered*) –

**shutdown** ()

Clean up and close connection to host computer. All control will be disabled. This can be called manually, but will also run automatically when your script exits.





**b**

`blue_interface`, 1



## B

`blue_interface` (*module*), 1  
`BlueInterface` (*class in blue\_interface*), 1

## C

`calibrate_gripper()` (*blue\_interface.BlueInterface method*), 1  
`cancel_gripper_command()` (*blue\_interface.BlueInterface method*), 1  
`command_gripper()` (*blue\_interface.BlueInterface method*), 1

## D

`disable_control()` (*blue\_interface.BlueInterface method*), 1  
`disable_gripper()` (*blue\_interface.BlueInterface method*), 1

## E

`enable_gripper()` (*blue\_interface.BlueInterface method*), 1

## G

`get_cartesian_pose()` (*blue\_interface.BlueInterface method*), 1  
`get_gripper_effort()` (*blue\_interface.BlueInterface method*), 2  
`get_gripper_position()` (*blue\_interface.BlueInterface method*), 2  
`get_joint_positions()` (*blue\_interface.BlueInterface method*), 2  
`get_joint_torques()` (*blue\_interface.BlueInterface method*), 2  
`get_joint_velocities()` (*blue\_interface.BlueInterface method*), 2  
`gripper_enabled()` (*blue\_interface.BlueInterface method*), 2

## I

`inverse_kinematics()` (*blue\_interface.BlueInterface method*), 2

## S

`set_joint_positions()` (*blue\_interface.BlueInterface method*), 2  
`set_joint_torques()` (*blue\_interface.BlueInterface method*), 3  
`shutdown()` (*blue\_interface.BlueInterface method*), 3